

Commode instruction leaflet

Reissue

If you reissue or are about to reissue this product, please thoroughly check all components for their safety.

This includes:

- **PAT test for electrical items**
- **Tightness of all nuts / bolts / castors**
- **Other screw in / bolt in / push in components.**
- **Also check all upholstery for security, splits, etc.**

If in any doubt, please do not issue or use, but immediately contact your supplier for service support.

When reissuing, a risk assessment may be required as to the suitability of the chair for a particular user.

Maintenance


Please make a safety check of the product at regular intervals or if you have any concerns.

Do not allow children or unauthorised persons to play with the equipment or use without proper supervision.

Tel 01942 243999 Fax: 01942 217 992
SimplyMed Limited
Unit B, New Venture Business Park
Leopold Street, Wigan, WN5 8DH
Email: info@simplymed.co.uk
Web: www.simplymed.co.uk



SimplyMed
shaping mobility



Code Description Weight limit

SMT007 Fixed Height Commode 154kg/25 stone

SMT008 Adjustable Height Commode 154kg/25 stone

SMT009 Detachable Arm Commode 154kg/25 stone

1. All Commodes

Please examine your commode for any visible damage before using. If you see any damage or suspect a fault, please contact your supplier.

In this case, please do not use your product as it may compromise your safety.

2. Adjustment of Legs (if applicable)

Remove the 'E' clip from the extension leg and extend all legs evenly to the desired height ensuring you do not exceed the last punched hole of the extension leg. (see fig. 1)

Once you have achieved the height you need, re-install the 'E' clip ensuring that it has passed through the extension leg and is fitting neatly and snugly around the commode leg. (see fig. 2)

Make sure that all legs are sitting evenly onto the floor; never extend legs in a manner that the commode is slanting. This will compromise the safety of the user.

3. Surrey Commodes only

The Surrey Commode has detachable arms; to remove these simply lift the arms out the locating tubes welded to the side of the commode frame. Always ensure that the arms are re-fitted neatly and evenly into the locating tubes before use.

Failure to do so will compromise the safety of the user.

4. Commode Buckets (All Commodes)

To remove the commode bucket, simply lift the toilet seat and fix the lid provided onto the bucket to avoid splashing and carefully raise the bucket and dispose of the contents into a proper sluice or WC.

Never use an open drain, sink or basin. This could cause a serious health risk.

Clean out the residue with a suitable disinfectant and refit to the commode after use.

5. Cleaning (All Commodes)

All commodes must be cleaned using a non-abrasive cleaner or can be pressure washed. However, always ensure that any chemicals you use with your pressure washing will not affect chrome plate, polymer coating or epoxy powder coating. After cleaning use a water dispersing agent such as WD40 on all joints and castors to reduce the chance of premature rust or seizure.