



**SimplyMed**

# REHABILITATION SUPPLIES & AIDS FOR DAILY LIVING



MOBILITY

TOILETING

BATHING

HOMECARE

BARIATRIC

REHABILITATION

[www.simplymed.co.uk](http://www.simplymed.co.uk)



We are proud to announce the release of our new SimplyMed Catalogue and trust you will find the products that you are looking for inside.

SimplyMed are one of the UK's leading privately owned, independent suppliers of mobility products and ADL equipment.

We supply products to hundreds of NHS departments, councils and distributors across the UK and have established a reputation for providing exceptional customer service and superior quality products on all fronts.

At SimplyMed we are passionate about delivering customer service which is intentional, consistent and delivers added value to our clients ensuring a professional and reliable service at all times. We are committed to training our staff to enable them to deliver an exceptional standard of service from point of first contact through to delivery of goods.

SimplyMed Ltd is a member of the British Healthcare Trades Association (BHTA) which means we have signed up to a Code of Practice giving our customers the assurance they are dealing with a reputable Company.

With this in mind we work closely with prescribers and users to develop products that are innovative, of high quality and offer outstanding value for money at all times.

We look forward to working with you and being of service going forward.

The SimplyMed Team



## MOBILITY PRODUCTS

### ADJUSTABLE HEIGHT WHEELED WALKING FRAME (STANDARD HOSPITAL SIZE)

A top quality, top value aluminium wheeled walking frame designed to meet the requirements of hospitals, nursing homes and everyday home use. Maximum user weight 159kg (25st). Maximum width 570mm. Maximum depth 540mm.

Code	Size	Handgrip Height
SMM001	Small	700 - 775mm 28" - 31"
SMM002	Medium	790 - 865mm 31" - 34"
SMM003	Large	885 - 960mm 35" - 38"
SMM00KW	Children's	585 - 660mm 23" - 26"

### ULTRA NARROW WHEELED FRAME

A top quality, top value aluminium wheeled walking frame designed to meet the requirements of hospitals, nursing homes and everyday home use. Maximum user weight 159kg (25st). Maximum width 490mm. Maximum depth 530mm.

Code	Size	Handgrip Height
SMM007	Small	700 - 775mm 28" - 31"
SMM008	Medium	790 - 865mm 31" - 34"
SMM009	Large	885 - 960mm 35" - 38"
SMM050	Aluminium Frame Stubs with Type Z Ferrule (pair)	
SMM051	Walking Frame Wheels (pair)	

- Moulded handles for improved infection control and comfort
- Sturdy yet soft wheels for added comfort
- No internal screw heads to prevent patient injury
- Double cross top bar to aid stability
- Anodized finish to improve aesthetics, cleanliness and longevity

### REPLACEMENT FERRULE

Code	Size
SMM055	25mm



AVAILABLE THROUGH APPROVED SUPPLY CHAIN ROUTE



## ADJUSTABLE HEIGHT UNWHEELED WALKING FRAME (STANDARD HOSPITAL SIZE)

A top quality aluminium unwheeled walking frame designed to meet the requirements of hospitals, nursing homes and everyday home use. Maximum user weight 159kg (25st). Maximum width 570mm. Maximum depth 540mm.

Code	Size	Handgrip Height
SMM004	Small	27" - 30" 685 - 760mm
SMM005	Medium	30" - 33" 760 - 835mm
SMM006	Large	34" - 37" 865 - 940mm
SMM00K	Children's	21" - 25" 535 - 635mm

## ULTRA NARROW UNWHEELED FRAME

A top quality aluminium unwheeled walking frame designed to meet the requirements of hospitals, nursing homes and everyday home use. Maximum user weight 159kg (25st). Maximum width 490mm. Maximum depth 530mm.

Code	Size	Handgrip Height
SMM010	Small	27" - 30" 685 - 760mm
SMM011	Medium	30" - 33" 760 - 835mm
SMM012	Large	34" - 37" 865 - 940mm

## RUBBER FERRULES

A range of high quality moulded rubber ferrules with metal inserts.

Code	
SMM052	Type Z 16mm
SMM053	Type Z 18mm
SMM054	Type Z 22mm
SMM055	Type Z 25mm
SMM052B	Type D 18mm
SMM053B	Type D 19mm
SMM054B	Type D 22mm



## GLIDE DISCS

Glide discs are attached to the rear ferrules to aid smooth movement on carpet.

Code	
SMGD002	Pair



## COMPACT TROUGH FRAME

Fully adjustable in height, handle length and handle angle, our Compact Trough Frame features ultra smooth 'PU' arm pads for added comfort and pressure sore prevention, as well as improved infection control. With optimum frame width this product can be used in hospital or home environment.

- Features robust steel troughs lined with ultra smooth PU for comfort and prevention of pressure sores
- Handles are adjustable for both length and angle
- Max user weight 160kg (25st)
- Tip Size 25mm
- Replacement pads available

Code	Size	Handgrip Height
SMM016	Small	810 - 885mm
SMM017	Medium	900 - 975mm
SMM018	Large	1000 - 1075mm
SMM018L	Extra Large	1200 - 1275mm
SMM122	Replacement Pads	



SMM122



## ELBOW CRUTCH - COMFY HANDLE

Comfy Handled Crutches have been designed to offer an even spread of pressure over the whole palm area offering more comfort for the longer term user over traditional crutches.

- Product Weight 715g
- Maximum User Weight 160kg (25st)
- 22mm ferrule (SMM054)
- Cuff dia 95mm



Code	Size	Handgrip Height	Cuff to handgrip
SMM028	Medium	23-33"	8"-11"
SMM029	Adult	27-37"	8"-11"



## ELBOW CRUTCH - PVC HANDLE

Used extensively in the NHS this versatile inexpensive and robust crutch. Adjustable between cuff and handle and handle and floor.

- Cuff diameter 95mm
- Maximum user weight 160kg (25st)



Code	Size	Handgrip Height	Cuff to handgrip
SMM019	Adult	25-36" (635-914mm)	9-12" (229-305mm)
SMM020	Medium	21-31" 533-787mm	8-11" 203-279mm
SMM021	Small	17-26" (432-660mm)	7-9" (179-229mm)
SMM022	X Long	32-42" (813-1067mm)	9-12" (229-305mm)



BEST SELLER



## ELBOW CRUTCH - PLASTIC HANDLE

The shape of the handle on this elbow crutch has been designed to spread the pressure more evenly over the hand, and the triangular shape and angle of the handle reflects the natural shape of the hand. Single adjustable crutches are adjustable between handle and ground. Double adjustable crutches are adjustable between cuff and handle and handle and ground. Latex free (including ferrule).

- Cuff Diameter 95mm
- 22mm Ferrule (SMM054)
- Maximum User Weight 180kg (28st)



Code	Description	Handgrip Height	Cuff to handgrip
SMM020P	Double adj - medium	21-31"	8-11"
SMM019P	Double adj - adult	25-36"	9-12"
SMM022P	Double adj - X long	32-42"	9-12"
SMM0225P	Single adj - adult	26-36"	9"



## CHILDREN'S ELBOW CRUTCH

High strength, lightweight double adjustable crutch with a soft touch, comfortable handle.

- Maximum User Weight 102kg (16st)
- 16mm ferrule (SMM052)
- Cuff dia 65mm (2½")

Code	Description	Handgrip Height
SMM023	4 - 7 years	380-530mm (15-21")
SMM024	6 - 10 years	440-670mm (17-26")



## AXILLA CRUTCHES

Traditional under arm crutch made from lightweight aluminium, simple to adjust and easy to clean

- Double adjustable
- Made from lightweight aluminium
- High density rubber
- Maximum User Weight 127kg (20st)
- 22mm ferrule (SMM054)



Code	Size	Handgrip Height
SMM025	Small	940-1145mm (37-45")
SMM026	Medium	1120-1320mm (44-52")
SMM027	Large	1320-1520mm (52-60")



## ADJUSTABLE FISCHER STICKS - LEFT & RIGHT HAND

A unique, biometrically shaped handle designed to evenly spread pressure across the palm of the hand. Useful for arthritics or those with poor grip. Available in various sizes.

- Maximum user weight 160kg (25st)

Code	Size	Height
SMM034	Medium Left	27"-36.5" (686-927mm)
SMM035	Medium Right	27"-36.5" (686-927mm)
SMM036	Large Left	32"-42" (813-1067mm)
SMM037	Large Right	32"-42" (813-1067mm)



## ALUMINIUM SUPPORT STICKS

A strong, lightweight support stick which is height adjustable.

- Maximum user weight 160kg (25st)

Code	Size	Height
SMM038	Medium	28" - 38" (711-965mm)
SMM039	Large	33" - 43" (838-1092mm)



**BEST SELLER**

MAXIMUM USER WEIGHT 160KG/25ST

## TRIPODS AND TETRAPODS

The centre of the handle on these sticks is barrel shaped to fit the palm. The front tapers to accommodate the fingers and encourage the user to point the index finger along the stem, and then flares at the end to prevent slipping.

- Maximum user weight 127kg (20st)

Code	Description
SMM044	Tetrapod Small Base
SMM045	Tetrapod Large Base
SMM046	Tripod
SMM047	Tetrapod CTR Legs



## LIGHTWEIGHT ROLLATOR WITH BAG

Light to carry, this New Aluminium rollator offers many upmarket features including PU seat for durability, 7" wheels for stability and a discreet carry bag.

### Features

- Very lightweight at 7kg or less
- Height adjustable handles 32-39"
- Arthritic-friendly loop cable brakes
- Underseat storage bag
- 7" / 18cm diameter solid tyres
- Seat Width 36cm (14")
- Seat Height 55.8cm (22")
- Weight Capacity 133kg (21st)

Code	Description
SMM049	Lightweight Rollator with Bag



**BEST SELLER**

## TRI WALKER WITH BAG

An adjustable height 3 wheel walker with easy to operate loop brakes. Utilises specially designed 8" molded wheels with extra width to avoid drain entrapment. Folds away easily for transportation and storage. Includes PVC shopping bag.

- Maximum user weight 133kg (21st)
- Height adjustment 33" - 38"
- Overall width 27"
- Overall depth 24"

Code	Description
SMM048	Tri Walker with Bag
SMM048B	Tri Walker Bag



SMM048B

## EXTRA TALL TRI WALKER

Our adjustable height Tri Walker as SM048 but with modified handles to give an extra 4" height on the handles giving a handle height of 37" to 42".

- Maximum user weight 133kg (21st)

Code	Description
SMM048T	Extra Tall Tri Walker



AVAILABLE THROUGH APPROVED SUPPLY CHAIN ROUTE



## SKANDIA, SINGLE & DOUBLE FLAP

The Skandia is a coated toilet frame which is designed to fit over the majority of WC suites. The Skandia offers additional height and integral arms to aid transfer. There are also the options of a floor fixing model or a separate kit for permanent installation if required. The unit is supplied complete with a splash guard.

- Easy height adjustment with horseshoe clips
- High quality toilet seat
- Removable splash guard
- Floor fixing kit available
- Integral arms to aid transfer
- Ideal for reissue
- Maximum user weight 127kg (20st)

Code	Description
SMT001	Skandia
SMT036	10LTR Bucket
SMT004	Sleeve for Skandia
SMT020	Skandia Double Flap (shown)
SMT016	Floor Fixing Kit Angled



**BEST  
SELLER**



## RAISED TOILET SEAT AND FRAME

A sturdy height adjustable, coated metal frame with a clip-on, one piece moulded seat. The seat has been designed with the comfort of the user in mind and features a partial cut-away at the front and rear to facilitate personal hygiene. The whole frame can then be positioned over the toilet.

- Width between arms 495mm
- Depth 480mm
- Adjustable height seat 450-600mm

Code	Description
SMT018	Raised Toilet Seat and Frame (35st)
SMT018W	21" Wide Raised Toilet Seat and Frame (41st)

## TOILET FRAME

Adjustable height, plastic coated toilet surround with rigid double strut steel frame.

Moulded plastic arm rests.

- Frame width 480mm (19")
- Depth 430mm (17")
- Height adjustment range 635-785mm (25" - 31")
- Maximum user weight 190kg (30st)
- Width adjustment (SMT014) 480-610mm (19" - 24")

Code	Description
SMT015	19" Toilet Surround
SMT014	Adjustable Height & Width Toilet Surround
SMT017	Floor Fixing Kit (not shown)



**BEST  
SELLER**

## ECONOMY COMMUNE WITH 5 LITRE BUCKET

Hygienic, easy-clean and comfortable commune chair has a fixed flame retardant padded back and detachable over seat cover with a high density polypropylene pan and lid. Splayed leg design provides optimum stability and allows them to be stacked for storage.

- Maximum user weight 159kgs (25st)
- Width between arms 400mm
- Depth of Seat 360mm
- Seat Height 480-570mm (adjustable version)

Code	Description
SMT007	Commune Standard
SMT008	Commune Height Adjustable
SMT009	Height adjustable stacking commune with detachable arms



**BEST  
SELLER**



DETACHABLE  
ARMS



FIXED  
HEIGHT

## MOBILE COMMUNE WITH FOOTRESTS

Chrome plated mobile commune which is able to glide over most raised toilet seats and clear most cisterns. Chrome plated finish and vinyl welded upholstery. The commune seat is made from blow moulded ABS plastic.

- Detachable arms
- Four brake castors
- Pram handle
- Colour: Grey, chrome
- Width between arms 450mm
- Depth of Seat 470mm
- Maximum user weight 165kgs (26st)

Code	Description
SMT006	Economy Mobile Commune
SMT028	Adjustable Height Mobile Commune
SMT006L	Replacement Footrests (pair)



SMT006L  
REPLACEMENT  
FOOTRESTS



## MALE URINAL

Male urinal with handle and anti-spill clip on lid

Code	Description
SMT049	1000ml Male Urinal



## FEMALE URINAL

Convenient and portable with carry handles and clip on lid

Code	Description
SMT048	1000ml Female Urinal



## 5 & 20 LITRE BUCKET

Replacement 5 & 20 litre bucket for use on most standard commodes

Code	Description
SMT003	5 Litre Bucket
SMT026	20 Litre Bucket



## ERGONOMICALLY DESIGNED ABS PLASTIC RAISED TOILET SEAT

Convenient, ergonomically designed raised toilet seat. It comes in either 50mm, 100mm or 150mm raise and is designed to fit the vast majority of UK toilet bowl shapes. The side mounted brackets have rubber over-moulds and the front skirt has a rubberised pad for maximum grip and stability.

- Ergonomically designed for comfort
- Smooth aperture interior to reduce waste build up
- Large aperture and front cutaway
- Maximum user weight 190kgs (30st)

Code	Description
SMT010	50mm (2") Raised Toilet Seat
SMT011	100mm (4") Raised Toilet Seat
SMT021	150mm (6") Raised Toilet Seat



## BATHBOARD/SHOWERBOARD

**BATHBOARD:** Made of highly durable reinforced PVC the Bathboard is compatible with plastic, steel and cast iron baths and is secured in position with a simple self locating system. 26", 27" & 28" widths available.

**SHOWERBOARD:** A strong showerboard made of highly durable reinforced PVC for maximum durability. Secured into position with a unique self-locating fixing system.

- Maximum user weight 190kg (30st)

Code	Description
SMB019	26" Bathboard
SMB020	27" Bathboard
SMB021	28" Bathboard
SMB022	26" Showerboard
SMB023	27" Showerboard
SMB024	28" Showerboard



## BATHSEAT

A strong bathseat made from reinforced PVC. Fixes securely with large suction pads and easy release tabs. Finished in brilliant white with soft grey trim. 6", 8" & 12" heights available. Meets MHRA Guidelines on entrapment.

- Maximum user weight 190kg (30st)

Code	Description
SMB025	6" Bathseat
SMB026	8" Bathseat
SMB027	12" Bathseat



## SHOWER CHAIRS

Tough plastic durable coating. Manufactured from strong tubular steel, featuring a stackable design. The comfortable polypropylene moulded seat has been designed with the safety of the user in mind, and has 5 drainage holes of 6mm diameter.

- Maximum user weight 127kg (20st)

Code	Description	Seat width	Height
SMB001	Fixed height	455mm	480mm
SMB002	Adj Height,	450mm	460-630mm
SMB003	Mobile	460mm	510mm
SMB003A	Adj Height, mobile	510mm	600mm



## HEIGHT ADJUSTABLE SHOWER STOOL

Our shower stool is a high quality shower stool designed for heavy duty use. The shower stool is adjustable in height and features sturdy, integral side handles to aid transfer. The shower stool has the added benefit of a clip-on/clip-off seat, making it simple to clean and reissue if required.

- Clip-on/ clip-off plastic seat
- Weight limit 159kg (25st)
- Height of seat from floor 495mm - 645mm
- Depth of Seat 275mm
- Width of Seat 370mm

Code	Description
SMB004	Height Adjustable Shower Stool



## SHOWER SEAT

This standard shower seat offers a competitively priced solution without loss of quality or function. The seat is height adjustable (in 25mm increments) enabling it to be set at various heights to suit the clients needs. The legs have fine height adjustments to compensate for any variations in the shower tray. The frame is manufactured from first grade 6,000 series alloy for greater longevity and the injection moulded seat can be clipped on and off for easy cleaning.

- **Maximum user weight 190kgs (30st)**
- **Wall fixed with adjustable leg support**
- **Simple to install, clean and maintain**

Code	Description
SMB006	Shower Seat
SMB007	Shower Seat With Padded Top



## DELUXE SHOWER SEAT

The Deluxe Shower Seat with Arms is height adjustable and the legs have fine adjustments to compensate for any variations in the shower tray. The integrated arms lift independently enabling safe and controlled transfer. The shower seat's frame is manufactured from rust free aluminium and the clip-on seat is injection moulded.

Code	Description
SMB008	Standard - No Padding
SMB009	With Padded Back and Seat



## CORNER SHOWER STOOL

An adjustable height corner shower stool which sits in the corner of the shower base. Polypropylene seat with drainage slats, removable for ease of cleaning. Warm touch, high performance plastic coating.

- **Maximum user weight 158kgs (25st)**
- **Height adjustment range 495mm - 647mm (19½" - 25½")**

Code	Description
SMB005	Corner Shower Stool
SMB015	Corner Shower Stool (with PU seat)



## PERCHING STOOLS

A range of practical, heavy duty stools manufactured in preferred high grade tensile steel, with a tough peel-resistant plastic coating. All models have front to back bracing struts for additional strength. The range incorporates new P.U. moulded components which include ergonomically designed robust seats which are non-slip in a water environment, fire retardant, hygienic and suitable for cleaning in industrial/ commercial washing machines. The seats, arms and component backs contain hygienilac which makes them resistant to MRSA and other microbiological contaminants.

- **Maximum user weight 159kgs (25st)**

Code	Description
SMH003	Perching Stool
SMH004	Perching Stool with Arms
SMH005	Perching Stool with Arms and Plain Back
SMH006	Perching Stool with Arms and Padded Back



## VINYL PERCHING STOOLS

We have a range of practical, heavy duty perching stools, manufactured in high grade steel with a tough plastic coating. All models have front to back bracing struts for additional strength. The seats 13½" x 11½" (340mm x 280mm) are padded and angled for improved posture and are PVC heat sealed conforming to the required fire regulations. Extension legs are chrome plated to BS1224 and are fitted with non-slip ferrules. White or Brown seats available.

- **Maximum user weight 159kg (25st)**

Code	Description
SMH028	Perching Stool
SMH029	Perching Stool with Arms
SMH030	Perching Stool with Arms and Plain Back
SMH031	Perching Stool with Arms and Padded Back





## BACKREST

A quality quick and simple way to help you sit up in bed. Angle adjustable and foldable for easy storage.

- Maximum user weight 108kgs (17st)

Code	Description
SMH010	Backrest



## ADJUSTABLE TROLLEY WALKER

This market leader is a very robust, all welded plastic coated steel design with angled sides to provide safe hand grips and install user confidence. The strong shelves are plastic injection moulded and each has 4 raised edges and the bottom shelf is set forward for ease of walking and dining. No tools necessary to detach trays.

- Maximum user weight 125kgs (19st)
- Shelf size 310mm x 405mm
- Overall width 455mm
- Overall depth 480mm

Code	Description
SMH001	Adjustable Trolley Walker
TC004	Replacement Tray



## HIGH BACK SITTING CHAIR

Brown vinyl, metal frame, winged head support and lumbar support. Completely adjustable (overall height, seat depth and back height). Generously padded and covered in vinyl for easy cleaning. Materials meet current fire safety standards.

- Maximum user weight 140kgs (22st)
- Overall height 1120mm - 1250mm
- Width 625mm
- Overall depth 670mm
- Seat height 430mm - 560mm
- Seat depth 450mm
- Arm height 675mm - 800mm

Code	Description
SMH009	Sitting Chair with Wings
SMH008	Sitting Chair No Wings (shown)



## BED LEVER

This bed lever is designed to assist transfer into and out of one side of the bed without affecting other users of the bed. The hand grip sits above the bed enabling a safe and simple transfer, it can be adjusted in height to suit various needs. The base frame is large enough to provide excellent stability and a traversing strap is available if required to further increase safety. The hand grip has been developed to MHRA guidelines. The Bed Transfer Aid can also be used on slatted beds when used in conjunction with the strap and buckle step.

- Height adjustable grip with safety rail
- Large base frame for increased safety
- Additional traversing strap available for greater stability
- High quality epoxy resin finish
- Easy to clean
- Suitable for industrial cleaning
- Conforms to MHRA guidelines
- Suitable for slatted beds when used with safety strap

### SPECIFICATIONS

Width of rail (exposed section) 325mm  
Width of base (section under bed) 605mm  
Height of rail (exposed section) 580mm-730mm  
Depth of base (section under bed) 700mm

Code	Description
SMG009	Bed Lever
SMG010	Bed Lever with Strap
SMG029	Spare Strap



## OVERBED TABLE

Our versatile table can be used around the house and is ideal for dining, reading and writing on.

- Attractive beech effect
- Easy clean and wipe down surface
- Fine angle adjustment of table top



Code	Description
SMH011C	Overbed Table No Castors
SMH013C	Overbed Table 2 Castors
SMH012C	Overbed Table 4 Castors

## MATTRESS ELEVATORS

The SimplyMed mattress elevator is ideal for raising or lowering of the mattress whenever required. It can also help users with respiratory problems achieve a more comfortable sleeping position.

- Maximum user weight 190kgs (30st)
- Linak Actuator
- Soft touch hand control
- Adjustable from 10° - 70°

Code	Description
SMH042	Mattress Elevators



## BARIATRIC FRAMES

- Heavy duty steel construction
- Chromium plated finish
- Comfortable moulded hand grips
- Non-slip rubber feet
- Height adjustable
- Maximum user weight 318kgs (50st)

Code	Description	Handgrip Height
SMF010S	Small	63-71cm
SMF010M	Medium	74-82cm
SMF010L	Large	84-92cm



## BARIATRIC CRUTCH

- Heavy duty steel construction with oversized cuff for comfort
- Chromium plated finish
- Comfortable moulded hand grips
- Non-slip rubber feet
- Height adjustable
- Maximum user weight 318kgs (50st)

Code	Description	Handgrip Height
SMF006	Bariatric Crutch	66 - 91cm




## BARIATRIC ROLLATOR

Our rolling walker has an extra wide design and robust construction for comfort and stability.

- Height adjustable handle - 820-870mm
- Padded seat
- 8" wheels
- Max user weight 227kg (36st)

Code	Description
SMM061	Bariatric Rollator



## BARIATRIC FOLDING FRAME

Our bariatric folding frame is designed with durability and stability in mind.

- Height adjustable: 810-990mm
- Max width: 740mm
- Max user weight 227kg (36st)

Code	Description
SMF010F	Bariatric Folding Frame
SMF010FW	Wheels for Folding Frame





## BARIATRIC RTS AND FRAME

Sturdy, extra wide and height adjustable frame with armrests.

- **Adjustable seat height** 410 - 610mm
- **Width** 690mm
- **Maximum user weight** 190kg (30st)

Code	Description
SMF008	Bariatric RTS and Frame



MAXIMUM  
USER WEIGHT  
190KG/30ST

## SIMPLYMED BARIATRIC HIGH BACK CHAIR

The SimplyMed Bariatric High Back Chair is specially designed for bariatric users. The standard finish is easy-to-service brown vinyl, other colours and fabrics are available subject to quantity. The padded arms extend and wrap around the front of the frame to aid transfer and enhance comfort. The bariatric chair is wider and deeper for comfortable bariatric use and incorporates fold away transfer wheels for easy manoeuvrability.

### SPECIFICATION:

- 260kg (41 stone)
- Height adjustable
- Ergonomic padded arms
- Seat wings as standard
- Transfer wheels as standard

### ADDITIONAL SPECIFICATION

- Colour: Brown & Black
- Weight: 4.90 kg
- Width: 610 mm
- Depth: 495 mm

Code	Description
SMF003	Bariatric High Backed Chair



## BARIATRIC COMMODE

Our Bariatric Commodes have the benefit of two fold up wheels for easily transporting the commode. The white top clips on and off for easy cleaning, it also has an integral aperture. Bariatric commodes feature flared arms as standard, padded back and overseat and a 5 litre bucket with lid.

- **Seat height** 450 - 580mm
- **Width between arms** 735mm
- **Maximum user weight** 248kg (39st)

Code	Description
SMF016	Bariatric Commode



## BARIATRIC PERCHING STOOL WITH PADDED BACK AND ARMS

Ideal for bathroom, shower or kitchen use. Providing extra wide seated support for larger / heavier users.

- **Width: 710mm**
- **Seat Width: 510mm**
- **Seat height: 460-610mm**
- **Max user weight 279kg (44st)**

Code	Description
SMF007	Bariatric Perching Stool



## MEDILAST 50M

A cost effective system of Progressive Resistive Exercise bands. Colour coded in five strengths to help develop exercise programs that are

- **Comfortable**
- **Progressive**
- **Versatile**

All in a system that utilises a simple colour coding to identify levels of resistance, aiding clinicians use and patient progression through structured treatment.

**Also available in a latex free version.**

Code	Description	Colour
SMTR204	Medilast 50m	Yellow
SMTR205	Medilast 50m	Red
SMTR206	Medilast 50m	Green
SMTR207	Medilast 50m	Blue
SMTR208	Medilast 50m	Black

SMTR209	Medilast 50m	Yellow LF
SMTR210	Medilast 50m	Red LF
SMTR211	Medilast 50m	Green LF
SMTR212	Medilast 50m	Blue LF
SMTR213	Medilast 50m	Black LF



## UNIBAND LATEX FREE BAND

An alternative low cost band giving excellent performance at value prices. All 50 metre length.

Code	Colour
MED60	Yellow
MED70	Red
MED80	Green
MED90	Blue
MED100	Black





## EPICAST® DUAL COMPRESSION EPICONDYLITIS CLASP

Targeted dual action compression for fast effective pain

Code	Size	
SMTR038	X-Small	20-23cm (8"-9¼")
SMTR039	Small	23-25cm (9¼"-10")
SMTR040	Medium	25-28cm (10"-11")
SMTR041	Large	28-31cm (11¼"-12½")
SMTR042	X-Large	31-34cm (12½"-13½")



**FITS BOTH  
LEFT AND  
RIGHT ARM**

**BEST  
SELLER**

## CLASSIC PLUS ACUPUNCTURE NEEDLES

Commonly used and adopted as the basic standard acupuncture needles in NHS hospitals and pain clinics.

### Special Features:

- Silicon free
- Pain free insertion
- Copper coiled loop-end handle
- Top surgical grade stainless steel

Code	Description
SMTR001	Classic+ Needle 0.20 x 13mm
SMTR002	Classic+ Needle 0.20 x 25mm
SMTR003	Classic+ Needle 0.22 x 13mm
SMTR004	Classic+ Needle 0.22 x 25mm
SMTR005	Classic+ Needle 0.22 x 40mm
SMTR006	Classic+ Needle 0.25 x 25mm
SMTR007	Classic+ Needle 0.25 x 40mm
SMTR008	Classic+ Needle 0.25 x 50mm
SMTR009	Classic+ Needle 0.25 x 70mm
SMTR011	Classic+ Needle 0.28 x 40mm
SMTR012	Classic+ Needle 0.30 x 25mm
SMTR013	Classic+ Needle 0.30 x 40mm
SMTR014	Classic+ Needle 0.30 x 50mm
SMTR015	Classic+ Needle 0.30 x 70mm

When ordering please state needle size



## DUAL CHANNEL TENS

### Technical Specifications

Channel:	Dual, isolated between channels
Pulse Intensity:	Adjustable 0-80mA peak into 500 Ω load each channel, constant current
Pulse Rate:	2Hz - 120Hz (adjustable).
Pulse Width:	40µs - 260µs (adjustable).

### Function Modes

- B: Cycle Bursts, 2Bursts/sec. 9Pulses/Burst. 100Hz, width is adjustable
- C: Continuous mode. Pulse rate, pulse width and intensity are adjustable
- M: Modulated Width. Pulse width is automatically varied in an interval of 6 seconds. The modulation range of pulse width is from setting value to 35% less than the control setting value, then returns to the setting value. Rate, width and intensity are fully adjustable.

Code	Description
SMTR114	Dual Channel TENS Kit
SMTR116	Replacement Cable (1)



## DUAL CHANNEL DIGITAL TENS

### Technical Specifications

Channel:	Dual, isolated between channels
Pulse Intensity:	Adjustable 0-80mA peak into 500 Ω load each channel, constant current
Pulse Rate:	1Hz - 160Hz (adjustable), 1Hz/step (1-20Hz), 5Hz/step (20-160Hz)
Pulse Width:	50µs - 260µs (adjustable).
Patient Compliance Meter:	Shows the treatment times in hours
Timer:	15, 30, 60 minute and continuous mode selectable

### Function Modes

- B: Cycle Bursts, 2Bursts/sec. 9Pulses/Burst. 100Hz, width is adjustable
- C: Continuous mode. Pulse rate, pulse width and intensity are adjustable
- M: Modulated Width. Pulse width is automatically varied in an interval of 6 seconds. The modulation range of pulse width is from setting value to 35% less than the control setting value, then returns to the setting value. Rate, width and intensity are fully adjustable.

Code	Description
SMTR115	Dual Channel Digital TENS Kit



## SIMPLYMED BUDGET ELECTRODES

Code	Description
SMTR016	1 pack 5 x 5cm - 4 pads per pack
	10 packs 5 x 5cm - 4 pads per pack
	50 packs 5 x 5cm - 4 pads per pack
	100 packs 5 x 5cm - 4 pads per pack
SMTR017	1 pack 9 x 5cm - 4 pads per pack
	10 packs 9 x 5cm - 4 pads per pack
	50 packs 9 x 5cm - 4 pads per pack
	100 packs 9 x 5cm - 4 pads per pack







Unit B, New Venture Business Park  
Leopold Street, Wigan, Lancashire WN5 8DH  
Tel: 01942 243 999 | Fax: 01942 217 992  
Email: [info@simplymed.co.uk](mailto:info@simplymed.co.uk) | [www.simplymed.co.uk](http://www.simplymed.co.uk)